



*ASME 2010 First Global Congress on
NanoEngineering for Medicine and Biology (NEMB2010)*

Development of Micro-Spectrometer for Neural Probe with Wireless Power Feed

Uhn Lee

Gacheon University of Medicine and Science, Incheon, Korea

Sang H. Choi

NASA Langley Research Center, Hampton, VA

Kyo D. Song

Norfolk State University, Norfolk, VA

Yeonjoon Park

National Institute of Aerospace, Hampton, VA

NASA SAA #15546 Collaboration Activity



Deep Brain Stimulation (DBS) and its Limits

Performance:

Jolt suppression voltage only


- diagnosis and search of anomaly required
- a single function

Power Source:


Implanted battery with tether line

- painful and cumbersome
- battery lasts 3~5 years
- power-line vulnerable to disconnection

DBS Implantation



Implanted DBS





New Approach: Probe-Pin Device

Performance:


- Diagnosis and search of anomaly required
- Multi-functions integrated
 - Jolt suppression voltage
 - Neuro-chemistry by micro-spectrometer
 - Neuro-electricity
 - Brain temperature
 - Brain pressure

Power Source:

- Wirelessly powered thru human tissue
 - Micro-coil with train of magnetic pulses
 - Rectenna array for microwave coupling

Logic Circuit for Monitoring, Control, and Data Acquisition

Master Logic dependent PPD



Development of Micro-Spectrometer

- Negative Zone Plate (NZP)
- Characterization of NZP
- Breadboard level of Micro-Spectrometer
- Smart Optics Materials for Micro-Spectrometer
- Smart Optics Materials Characterization

Photon Intensity Simulation for Micro Ring Grating Spectrometer



❖ 3-D optical simulation for 500 μm ring grating (25 transparent rings) with Fresnel diffraction.



❖ Photon intensity map of optical axis (z-axis) vs. radius direction (x-axis) shows strong focusing effect and spectral separation.

❖ Spectral resolution ($\Delta\lambda$) is improved with the increasing number (n) of rings as $1/n$.

- $\Delta\lambda=20$ nm for 50 rings
- $\Delta\lambda=10$ nm for 100 rings

❖ Each wavelength has a strong focal point which is surrounded by two weak light cones.


❖ Circular aperture slit of 1% of the ring grating diameter will pass the photons of a selected wavelength.




Spectral scans:

Circular Grating: 100 rings, 757 μ m diameter
Aperture: 10 μ m diameter


Green Laser: 532nm



Red Laser: 633nm



Green & Red Lasers: 532nm & 633nm



Micro-Spectrometer Component Test Platform

Micro-Spectrometer

Pin Probe Device

Motor/Position Sensor Controller

Micro Linear Motor

Micro Position Sensor

Mini USB


XBEE
(ADC, Wireless Data Transfer)

InGaAs Detector



LEDs

OP-Amp



Switches



Configuration Diagram



Computer




Intensity vs. Z (distance)

$$\lambda = \frac{R_n^2}{n \cdot Z}$$


Intensity vs. λ (Wavelength)


[SPECTROMETER]

Micro-Spectrometer Component Test Platform Parameters




- Light Sources: White, Blue, and UV LED
- Integrated Light Pipe Design: Diameter = 3mm
- Micro Linear Motor: Travel Distance = 6mm
- Magnetic Micro-Position Sensor: 2 micrometer resolution
- Photo Detector: InGaAs Photo Diode
- Mini OP-Amplifier: Gain (V/I) = 10^7
- USB/Serial Interface to control linear motor
- Power: 3.3V, 1.8A max
- ADC & Wireless Data Transfer: XBEE Chipset (ZigBee Network)





μ -Spectrometer Development Plan





Smart Optics Materials




- Smart Optics Materials (SOM) are new optical materials that can control deep properties (intensity, phase, polarization, coherence) of passing light, including
 - electro-optic materials (non-linear optical crystals, liquid crystals, electro-optic polymers)
 - magneto-optic materials (Faraday Effect and Kerr Effect),
 - electro/thermo-chromic materials
 - chemicals that induce refractive index changes
 - optical materials that depend on temperature and pressure
- New prototype SOM characterization system with advanced software is under development at NASA LaRC.
 - Intensity, phase, and polarization of passing lights through many new optical materials while applying various physical/chemical quantities on the materials.
 - A miniaturized compact system has USB interface and exchangeable components for various R&D and commercial applications.
 - This innovative versatile SOM characterization system and characterized optical materials will be used for complete micro spectrometer system.
 - Commercialization of this system accelerates development of new optical materials and devices.

(a) Liquid crystal




Smart Optics Material Characterization System


(a) Electro-Optic Materials: liquid crystal, non-linear crystal




(b) Reflective samples: MEMS device, piezoelectric actuator materials, thermal expansion coefficient measurement



(d) Chemical sample holder, temperature/pressure control





(c) Magneto-optical materials: Kerr effect, Faraday effect











Smart Optics Material Characterization System





SOM Characterizer Software





Progress: | Cursor (418,241) | Offset (469,289) [51,48] | Mouse (693,26) | Ready



Dynamic Phase Map Analysis





Liquid Crystal Single Layer



Liquid Crystal Double Layer


Phase Angle Measurement






Mobile Wireless Power Transfer with Micro Coil and Microwave

Micro Induction Coil (MIC)




Flexible Membrane Inductor



Thin-Film Rectenna Array (TFRA)



3 x 3 flexible rectenna




Short Range (1~10cm) Solution




Wireless Power Transfer

Inductance Power Transfer


- 0-3 cm, short range
- Safe for human interaction
- low power applications




Planar Coils



Coil powering an LED



Coil Test Setup



8 mm TDS 2022B - 3:01:24 PM 8/25/2008

Microwave Power Transfer

- 1-1000 m, long range
- Low to high power applications






Flexible Rectenna Array




Microwave Generator and Amplifier

Thin-Film Rectenna Array (TFRA)




PPD Logic (internal) with TFRA or MIC





Features of Probe Pin Device

- Selection of Targeted Spectral Line for μ -SM:
 $600 \text{ nm} < \lambda < 1000 \text{ nm}$
- Spectral Resolution: $\leq 10 \text{ nm}$
- Downsize of Micro-Spectrometer:
 $< 2 \text{ mm}$
- Platform Design for Other Sensors Integration:
 ρ , T, V, and chemicals
- Wireless Power Feeder and Integration:
both MIC and MW
- Telemetry Circuit and Test with Skin Tissue:
FM mode
- Chip-scale Design of Electronics:
Miniaturization





Summary

- Neural Sensing Probe Pin Device with micro-spectrometer is under development.
 - Micro Spectrometer Component Test Platform was set up.
 - The world's smallest optical path volume (1mm^3 was achieved with 0.75mm diameter.)
 - This is $1/1000^{\text{th}}$ volume compared with today's smallest commercial spectrometers (1cm^3)
 - Demo model developed with actuator-controlled focal length as a stop-gap approach
 - All basic electronics for sensor, telemetry, and power coupling were fabricated, but required chip-scale miniaturization
- Prototype SOM characterization system is developed to measure **intensity, phase, and polarization** of passing lights through many new optical materials while applying various physical/chemical quantities on the materials.
 - A miniaturized compact system has **USB interface** and **exchangeable components** for various R&D and commercial applications.
 - SOM system and software are **Tech-Transfer ready**
- Wireless Power Transmission test using polyurethane layers and pig skin was satisfactorily performed.
 - MIC was proven to be effective for short range power transmission, like a hat system
 - TFRA was effective for short (near-field capacitive coupling) and long (far-field) range power transmission.